



Optimize Your Research

Recombinant Human cGMP-Dependent Protein Kinase 1 , C-6His, Human cells protein

Code CD02047

Storage: This lyophilized preparation is stable at 2-8°C, but should be kept at -20°C for long term storage, preferably desiccated. Upon reconstitution, the preparation is stable for up to one week at 2-8°C. For maximal stability, apportion the reconstituted preparation into working aliquots and store at -20°C to -70°C. Avoid repeated freeze/thaw cycles.

Intended Use

This product is for research use only, not for use in diagnosis procedures. It is highly recommended to read this instruction entirely before the use.

Source

Human Cells

Purity

Greater than 80% as determined by reducing SDS-PAGE.

Physical Appearance

Sterile Filtered White lyophilized (freeze-dried) powder.

Formulation

Supplied as a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, 0.25mM DTT, 0.1mM EDTA, 0.1mM PMSF, 20% Glycerol, pH 8.0.

Endotoxin

Less than 0.1 ng/µg (1 EU/µg) as determined by LAL test.

If you have any question on order please contact us via: order@bt-laboratory.com; technical assistance please contact us via: support@bt-laboratory.com More product visit www.bt-laboratory.com

5/F 2 Bldg., 501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China
Tel: 86 21 31007137 | E-mail: support@bt-laboratory.com | www.bt-laboratory.com