



Optimize Your Research

Recombinant Bovine Enterokinase Light Chain protein

Code CD02426

Storage: This lyophilized preparation is stable at 2-8°C, but should be kept at -20°C for long term storage, preferably desiccated. Upon reconstitution, the preparation is stable for up to one week at 2-8°C. For maximal stability, apportion the reconstituted preparation into working aliquots and store at -20°C to -70°C. Avoid repeated freeze/thaw cycles.

Intended Use

This product is for research use only, not for use in diagnosis procedures. It is highly recommended to read this instruction entirely before the use.

Source

Escherichia coli.

Molecular Weight

Approximately 28 kDa, a single non-glycosylated polypeptide chain containing 235 amino acids.

Purity

Physical Appearance

Sterile Filtered White lyophilized (freeze-dried) powder.

Formulation

50 mM Tris-HCl, pH 8.0, 0.5 M NaCl and 50 % glycerol.

Endotoxin

Less than 1 EU/μg of rBoEKL as determined by LAL method.

If you have any question on order please contact us via: order@bt-laboratory.com; technical assistance please contact us via: support@bt-laboratory.com More product visit www.bt-laboratory.com

5/F 2 Bldg., 501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China
Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com